

Preparing YPD media

Contents

Title : Liquid YPD media	1
Date	1
Results	1
Conclusions	1

List of Figures

Title : Liquid YPD media

Date

04-03-2021 ## Method

- Final Volume: 200 ml liquid YPD

Amount	Type
180 ml	H2O
2 g	Yeast extract
4 g	bacto-peptone
20 ml	20% dextrose

- Filter sterilize before adding dextrose
- Preparing 200ml of 20% dextrose

Amount	Type
200 ml	H2O
40 g	dextrose

Results

- how was the media checked?

Conclusions