

# Preparing YPD media

## Contents

<b>Title : Liquid YPD media</b>	<b>1</b>
Date . . . . .	1
Method . . . . .	1
Results . . . . .	1
Conclusions . . . . .	2

## List of Figures

### Title : Liquid YPD media

#### Date

04-03-2021

#### Method

- Final Volume: 200 ml liquid YPD

Amount	Type
180 ml	H2O
2 g	Yeast extract
4 g	bacto-peptone
20 ml	20% dextrose

- Filter sterilize before adding dextrose
- Preparing 200ml of 20% dextrose

Amount	Type
200 ml	H2O
40 g	dextrose

#### Results

- how was the media checked?

## Conclusions